

9.3 Replaceable Spare Parts List

When ordering spare parts, please be certain to supply the model number and serial number of your analyzer.

PART NO.	PART DESCRIPTION
	<u>Printed Circuit Boards</u>
10424310	Oxygen Board
10424570	Moisture Board
54000001	Front Display
	<u>Hardware Items</u>
SNT0010	Oxygen Sensor
63000308	Vacuum Pump – 110 VAC
63000309	Vacuum Pump – 240VAC
63000004	Vacuum Pump Rebuild Kit
76000018	Vacuum Pump Cleaning Fluid
65000100	Vacuum control needle valve
16233870	Purifier
15031560	Gas Panel Ass'y with heater – w/o purifier
16016910	Cap, oxygen sensor
60300268	VCR Filter Gasket
60300241	VCR Gasket
45002363	Fuse 6.3 Amps (110VAC)
45002331	Fuse 3.15 Amps (220VAC)
59017237	Power Cord – 110 VAC

Table 12: Spare Parts List